



Abilia

Lightwriter SL50




ABILIA

Lightwriter SL50

- Easy to use, fast and **dedicated keyboard-based (text-to-speech)** device for literate users
- Communicate efficiently either through face-to-face communication, SMS texting, or mobile phone calls
- It can be used both by typing the keys and by scanning with 1 or 2 control switches




Highlights of Value and Functionality

 Fast start-up and simple to use menu interface

- Personalized vocabulary for quicker communication


 Connect to a mobile phone using Bluetooth

 High quality Acapela voices or My Own Voice input

 Simple 'plug and play' 1-switch or 2-switch scanning

- Dual displays and outward facing speaker

 Backlit screens and keys which automatically adjust to different light levels

 Remote software updates



What is new with Lightwriter SL50?




- Fast start-up



- User 'Structured' Menus



- Backlit colour display



- Plug and Play Scanning

- Backlit keyboard



- Fast charge (4 hours) and all day battery life



- Weight 740g (>100g lighter)



- Remote SW updates via Ethernet port



- A selection of Acapela voices and import of new voices incl. voice banking



- Multi-lingual



- Import of personalized vocabularies



- Mobile phone and SMS (via Bluetooth)



Lightwriter SL50 vs Allora 2

USPs Lightwriter SL50

- Lighter (740 gr vs 1 kg or 1.6 lbs vs 2.2 lbs)
- Smaller
- Remote software updates via ethernet port
- Backlit screens and keyboard – can be operated in the dark
- 3 modes of scanning (Allora 2 modes)
- Connect to a mobile phone using Bluetooth




Lightwriter SL50



Allora 2

Lightwriter SL50 vs software app's

- **Lightwriter** is a dedicated product and can be used all the time; **“always ready to use”**
- The user gets help to guide his/her fingers to the correct keys by **keyguard**
- **Tactile feedback** via the keys
- **Water resistant**
- **Impact resistant**



Lightwriter vs tablet based AAC

- **The SL50 compared to a tablet is:**
 - More robust – designed for purpose
 - Less sensitive to drool and dribble
 - Easier to use for persons with tremors/shaking
 - Quick to get started with and easy to learn - turn on and type
- Not all users need or want a tablet-based device:
 - **Elderly:** some not used to tablet technology
 - **MND/ALS & Multiple Systems Atrophy:** user might struggle with tablet when hands become weak. SL50 keyguard gives better support and user can rest hands on keyguard without touching buttons
 - **Parkinson & Huntington:** SL50 has built in sensitivity controls to prevent double press caused by tremors. The keyguard give support for hands
 - **ASD/Autism:** SL50 is a dedicated device that can only be used for communication, a tablet filled with other applications can be distracting

