

# Turnplate, Turn Cushion, One-Way Slide



## NURSING AIDS

TURNPLATE SEAT-TURN CUSHION ONE-WAY SLIDE CUSHION BED ONE-WAY SLIDE



TURNPLATES For controlled rotation

Samarit Turnplates have already proven themselves tens of thousand times. They are featured by the handiness, controlls rotation, ease of cleaning, safety and longevity. Application in clinics, nursing homes and homecare.



#### TURNPLATE Blue 30

Size: Ø 30 cm / 40 cm Ø 12 inch / 16 inch ABS/Teflon, blue

Order No. 480.0100 UDI/GTIN 7640144461089

### SEAT & CHAIR TURN CUSHION ONE-WAY SLIDE CUSHION

For standing and seated transfer



#### SEAT & CHAIR TURN-CUSHION Blue 38

Size: Ø 38 cm / Ø 15 inch Antimicrobic surface, blue Biocompatibility DIN 10993-5 Flameretardant

Order No. 480.1038 UDI/GTIN 7640144461102



TURNPLATE Blue 38

Size: Ø 38 / 49 cm Ø 15 inch / 19 inch ABS/Teflon, blue

Order No. 480.0380 UDI/GTIN 7640144461072



TURNPLATE Bariatric

Size: Ø 41 cm / 54 cm Ø 15 inch / 21 inch ABS/Teflon, blue

Order No. 480.0400 UDI/GTIN 7640144461096

The Samarit Seat-Turn Cushion, is also manufactured from textile materials an has a solid quilted rotating cushion of flame resistant material. The bottom portion is equipped to be slide resistant. Application on car seat, settee, chair, etc.

Two soft discs turnable attached with turning device, tissue and foam used from fire resistent materials, bottom disc with non-slip, washable 40°.

Raise-Up Cushion slides to one side only. Permits to reposition patient in chair with little effort.



#### SEAT ONE-WAY SLIDE CUSHION Blue

Size: 45 x 42 cm / 18 x 17 inch Tissue with one-slide, Top PU-textile grey coating

Order No. 450.4060 UDI/GTIN 7640144461034

#### **High quality standards** SAMARIT transfer aids have:

- CE mark
- Design patents
- 2 years warranty

CE

## www.samarit.com

Download your language from our website or ask your local distributor.



Samarit Medical AG Gewerbestrasse 12 CH-8132 Egg, Schweiz Tel.: +41 44 918 10 11 Fax: +41 44 918 24 39 e-mail: info@samarit.com - www.samarit.com